



FW-03



Blood and Infusion Warmer

BLOOD AND INFUUSION WARMER





the liquid which will be infusioned into human body by warming the infusion pipe, with high efficiency and easy use; dry groove warmer structure, adopt conventional infusion pipe, no need special consumables.





PROFESSIONAL MEDICAL INSTRUMENT

WE NELL WILL NEVER CHANCE UNTIL DESTROYMENTS TO YOUR EACH REACTIVATE WELL MILE NEVER CHANCE UNTIL DESTROY
OFFICE YOUR EACH REACTIVATE WILL WILL

PROTECT YOUR CACHHEALTHY NEVER CHANGE UNTIL DEATH



Temperature controlling technology

Adopts intellectual microcomputer Controlling technology, with triple protection.

Infrared detection function can work continuously for 3 to 6 hours.








According to the different length and temperature requirements choose different ways of winding





3 kinds of choice











Blood and Infusion Warmer





| Work mode | Continuous work | Model | FW-03 |
|-----------|-----------------|------------|---------------------------|
| Power | ≤70W | Voltage | AC 220V ± 100% 50Hz ± 1Hz |
| Baro | 86kPa ~ 106kPa | Net weight | 0.8kg |

Packing List

| Power cable | X1 |
|-------------|----|
| Belt | X1 |
| Tube clamp | X2 |
| Manual | X1 |



Riransa Inc. Healthcare Limited

INDIA MARKETING OFFICE:

Head Office: 203, DDA Bulding No-2, District Centre, Janakpuri, New Delhi-58 **Corp. Office:** 61- Sharda Apartment, Ravindra Garden, Aliganj, Lucknow-226024

Mobile No: 8953961111

Website: www.riransa.net E-mail Id: info.riransa@gmail.com, care@riransa.net